




**EBERLE**

**Small Joint**



Control Unit







### Control unit

- Panel to control motor-powered EBERLE C3 handpieces
- 5-inch touch display
- Automatic brightness adjustment
- Softkeys for basic functions
- Adjustable oscillation
- Oscillation mode: revolutions or frequency
- User-specific assignment of the control elements for the footswitch as well as for the handpiece
- Individual acceleration behaviour
- User name editable
- Ergonomic 5-button footswitch

### Shaver System

- **Control unit**  
article no. 913055
- **Handpiece**  
article no. 913065SJ
- **Footswitch**  
article no. 913057

### Footswitch

article no. 913057

- user-friendly, ergonomic 5-button footswitch





### Small Joint Handpiece

article no. 913065SJ

- ergonomically shaped handpiece, autoclavable
- for operation of Small joint shaverblades SJ
- micro-processor controlled drive motor
- weight 185 g + 4 m cable
- speed ranges:
 

rotation	500 – 6000 min <sup>-1</sup>
oscillation	500 – 4000 min <sup>-1</sup>

### Cleaning Adapter

for Shaver-handpieces  
article no. 100020



for Universal Handpiece  
article no. 100025



for modular attachments  
article no. 100026




### Cleaning Brushes

for handpieces

**Cleaning brush, bent-over tip**  
length: 150 mm  
article no. 011001-301



**Cleaning brush with eyelet**  
length: 150 mm  
article no. 011001-302



**Cleaning brush, rounded tip**  
length: 300 mm  
article no. 011001-300



## Universal Handpiece manual operation

- ergonomically shaped handpiece, autoclavable
- micro-processor controlled drive motor
- manual operation of anticlockwise/clockwise rotation
- speed range: 0 – 1500 min<sup>-1</sup>, infinitely variable
- weight 750 g + 4 m cable



○ Universal Handpiece  
article no. 913075

..... AO Drill Adapter  
article no. 913044

..... Drill chuck adapter  
article no. 913027  
article no. 913019  
article no. 913022

..... Sagittal saw adapter  
article no. 913024

..... Wire Driver Adapter  
article no. 913026  
article no. 913034

### AO Drill Adapter

receptacle for standard AO  
article no. 913044



### Drill Chuck Adapter

receptacle for 0.7–7.4 mm  
article no. 913027



### Drill Chuck Adapter High Power

receptacle for 0.7–77.4 mm  
rotational speed 0–7500 min<sup>-1</sup>  
article no. 913019



### Drill Chuck Adapter

receptacle for 0–73 mm, no spanner  
article no. 913022



### Sagittal Saw Adapter

Oszillationsfrequenz 0–722.500 Zykl./min  
article no. 913024



### Wire Driver Adapter

receptacle for 1.8–3.2 mm  
article no. 913026

receptacle for 0.6–1.8 mm  
article no. 913034






## Shaver Blades for EBERLE Shaver System

for small joint surgery, effective length 70 mm  
reusable




Packing unit 5 pieces

	Ø 2.0 mm article no.	Ø 2.5 mm article no.	Ø 3.5 mm article no.
 <b>Full Radius Resector</b>	9131120SJ	9131125SJ	9131135SJ
 <b>Aggressive Full Radius Resector</b>	9131320SJ	9131325SJ	9131335SJ
 <b>Round Burr</b>		9131025SJ	9131035SJ
		Ø burr 1.8 mm	Ø burr 2.5 mm

### Cleaning Adapter

for reusable  
shaver blades  
article no. 100001



### Cleaning Set

article no. 900019






## Shaver Blades for Eberle Shaver System

for small joint surgery, effective length 70 mm  
sterile packaging, single-use



Packing unit 5 pieces

	Ø 2.0 mm article no.	Ø 2.5 mm article no.	Ø 3.5 mm article no.
 <b>Full Radius Resector</b>	9141120SJ	9141125SJ	9141135SJ
 <b>Aggressive Full Radius Resector</b>	9141320SJ	9141325SJ	9141335SJ
 <b>Round Burr</b>		9141025SJ	9141035SJ
		Ø burr 1.8 mm	Ø burr 2.5 mm

## Hook Probe, Small Joint

tip 3 mm, effective length 70 mm  
article no. 900355 <sup>1)</sup>





## Sagittal Saw Blades

Distributed by Eberle GmbH & Co. KG



	article no.
saw blade 25 x 6 x 0.6 mm cutting thickness 0.6 mm	900092



	article no.
saw blade 21 x 9 x 0.4 mm cutting thickness 0.4 mm	900116



	article no.
saw blade 31 x 11 x 0.4 cutting thickness 0.4 mm	900117



	article no.
saw blade 50 x 14 x 0.6 cutting thickness 0.6 mm	900119

## Wires

Distributed by Eberle GmbH & Co. KG



Packing unit 10 pieces

		article no.
Kirschnerdraht Trocar – round	Ø 0.8 x 150 mm	900108
Kirschnerdraht Trocar – round	Ø 1.0 x 150 mm	900107
Kirschnerdraht Trocar – round	Ø 1.2 x 150 mm	900106
Kirschner wire Trocar – round	Ø 1.2 x 160 mm	900082
Kirschner wire Trocar – round	Ø 1.2 x 310 mm	900112
Kirschner wire Trocar – round	Ø 1.4 x 150 mm	900105
Kirschner wire Trocar – round	Ø 1.6 x 150 mm	900090
Kirschner wire Trocar – round	Ø 1.8 x 150 mm	900089
Kirschner wire Trocar – round	Ø 2.0 x 150 mm	900088
Kirschner wire Trocar – round	Ø 2.0 x 310 mm	900111
Kirschner wire Trocar – round	Ø 2.4 x 310 mm	900086
Kirschner wire Trokar – Trokar	Ø 1.2 x 160 mm	900094
Kirschner wire Trokar – Trokar	Ø 1.2 x 310 mm	900084
Kirschner wire threaded trocar – trocar	Ø 1.2 x 200 mm	900095
Kirschner wire threaded trocar – round	Ø 1.2 x 200 mm	900083
Kirschner wire round – round	Ø 1.5 x 380 mm	900109
Kirschner wire round – round	Ø 1.6 x 230 mm	900028
Kirschner wire round – round	Ø 1.6 x 310 mm	900110

Packing unit 1 piece

		article no.
Drilling wire trocar, knurled with special eyelet	Ø 2.4 x 430 mm	900027
Drilling wire without eyelet	Ø 2.4 x 310 mm	900025
Drilling wire with special eyelet	Ø 2.4 x 380 mm	900026
Switching stick	Ø 4.0 x 230 mm	900085

## Endoscopes

autoclavable




	article no.
Arthroscope Ø 2.75 x 110 mm, 30°	900113 <sup>1)</sup>
Arthroscope Ø 2.75 x 110 mm, 30°, repair part replacement	900113R <sup>1)</sup>
Arthroscope Ø 1.9 x 60 mm, 30°	900309 <sup>1)</sup>
Arthroscope Ø 1.9 x 60 mm, 30°, repair part replacement	900309R <sup>1)</sup>
Arthroscope sterilization case, 290 x 80 x 52 mm	900091 <sup>1)</sup>

## Arthroscopy Shafts



	article no.
Arthroscope sheath with two rotatable taps, Ø 2.75 mm, 30°	900114
Arthroscope sheath with two rotatable taps, Ø 1.9 mm, 30°	900310

## Trocar, obturator




	article no.
Trocar pointed for article no. 900114	900115
Obturator conic blunt for article no. 900114	900120
Obturator conic blunt for article no. 900310	900311

<sup>1)</sup> distributed by Eberle GmbH & Co. KG

## Hand instruments

for small joint surgery




Distributed by Eberle GmbH & Co. KG




Tip profile

Bite width

Tip width

	tip profile	bite width	tip width	article no.
 <b>Basket punch,</b> junior, blunt, straight	1.70 mm	1.70 mm	3.00 mm	900326
 <b>Basket punch,</b> junior, oval, straight	1.20 mm	2.30 mm	3.50 mm	900327
 <b>Arthroscopy punch,</b> junior, straight	1.20 mm	2.30 mm	3.50 mm	900328



**WAS MACHEN WIR MIT DIESER FREIEN SEITE?**





**EBERLE**

*Made in Germany*

**EBERLE GmbH & Co. KG**

Dachsteinstraße 2  
75449 Wurmberg  
Germany

Phone +49 (0) 7044 9611-288

[sales@eberle-med.com](mailto:sales@eberle-med.com)  
[www.eberle-med.com](http://www.eberle-med.com)

