



EBERLE

HF


HF generator

distributed by Eberle GmbH & Co. KG

Compact HF generator for monopolar and bipolar arthroscopy and general surgery.
More indication profiles on request.



Description	Article no.
Electrosurgical unit EMED ATOM	900384



Description	Article no.
Footswitch for Electro-surgical unit EMED	900385
Footswitch for Electro-surgical unit EMED, wireless	900385-1

Examples for the application area



Knee joint

Partial meniscus resection, partial synovectomy, Plica synovialis resection, lateral release



Hip joint

Synovectomy, cartilage smoothing, labrum, partial resection



Shoulder joint

Acromioplasty, partial synovectomy, cartilage smoothing, tenotomy

Bipolar connection cables

distributed by Eberle GmbH & Co. KG



reusable

	Description	Article no.	Length
	Bipolar forceps cable, EMED SDS	900335	5 m
	Bipolar forceps cable, International 2-pole	900207	5 m
	Bipolar forceps cable, Martin Maxium	900374	5 m
	Bipolar forceps cable, Erbe VIO/ICC/ACC, Storz	900338	5 m

Other cables on request

Bipolar arthroscopy electrodes




distributed by Eberle GmbH & Co. KG

Safe handling, maximum efficiency.

The BiCart electrodes are used for cutting, vaporizing, coagulation and thermal modification of soft tissue.



reusable

		Shaft length = 115 mm	Shaft length = 170 mm
	Description	Bipolar vaporization electrode, "VAP"	Bipolar vaporization electrode, "VAP"
	Article no.	900197	900196
	Description	Bipolar vaporisation electrode 70°	Bipolar vaporisation electrode 70°
	Article no.	900195	900194
	Description	Bipolar hook electrode 90°	Bipolar hook electrode 90°
	Article no.	900193	900268
	Hook length	1.5 mm x Ø 0.8 mm	1.5 mm x Ø 0.85 mm

Bipolar forceps

distributed by Eberle GmbH & Co. KG



reusable

Description	Article no.
Bipolar HF forceps, tip 1 x 8 mm, 160 mm	900448
Bipolar HF forceps, tip 2 x 8 mm, 160 mm	900367

Monopolar HF handles, Evoline

distributed by Eberle GmbH & Co. KG

- For electrodes with 4 mm shaft
- Hexagonal anti-twist protection



reusable




HF handles, Evoline	U _{max}
HF handles for 4 mm shaft electrodes with integrated cables, two buttons	4.3 kVp



	Description	Monopolar HF handle, Erbe ICC/ACC, Storz
	Article no.	900419
	Length	3 m
	Description	Monopolar HF handle, Martin, Berchtold, Integra
	Article no.	900420
	Length	3 m
	Description	Monopolar HF handle, Valleylab, International 3 Pol
	Article no.	900421
	Length	3 m
	Description	Monopolar HF handle, EMED SDS
	Article no.	900435
	Length	3 m

Monopolar electrodes

distributed by Eberle GmbH & Co. KG

	Description	Article no.
	Neutral electrode cable EMED, 5 m	900410
	Disposable neutral electrode, self-adhesive	900395
	Silicone-Neutral electrode for adults, reusable, 4 m cable	900447





Active monopolar electrodes

distributed by Eberle GmbH & Co. KG

- Suitable for HF and RF generators
- Autoclavable
- Hexagonal anti-twist protection
- U_{max}: 4,3 kVp

 Ø 4 mm

reusable

	VPE = 5	
	Description	Knife electrode
	Article no.	900411
	Dimensions	26 mm 3.2 mm
	Description	Spatula electrode
	Article no.	900412
	Dimensions	22.5 mm 2.8 mm
	Description	Lancet electrode, strong blade
	Article no.	900413
	Dimensions	25 mm 2.5 mm
	Description	Knife electrode, lateral tip
	Article no.	900414
	Dimensions	24 mm 3.3 mm

Monopolar arthroscopy electrodes







distributed by Eberle GmbH & Co. KG

- Suitable for HF and RF generators
- Autoclavable
- Hexagonal anti-twist protection



Ø 4 mm

reusable

	Description	Article no.
	Arthroscopy electrode needle 90° angled, 1.5 x Ø 0.8 mm, 130 mm	900142
	Arthroscopy electrode needle 90° angled, 4.0 x Ø 0.8 mm x 130 mm	900143
	Arthroscopy electrode knife 45°, cranked, 4.0 x Ø 1.5 mm x 130 mm	900144
	Arthroscopy electrode pin 90° angled, 3.0 x Ø 1.5 mm, 130 mm	900145
	Arthroscopy electrode ball Ø 4.0 mm, 122 mm	900146
	Surface electrode Ø 3.5 mm, 70°, 115 mm	900442

Electrodes for plastic surgery




distributed by Eberle GmbH & Co. KG

- Suitable for HF and RF generators
- Autoclavable
- Hexagonal anti-twist protection
- Umax: 4,3 kVp



Ø 4 mm

reusable

	VPE = 5	
	Description	Ball electrode
	Article no.	900415
	Dimensions	Ø 1 mm
	Description	Ball electrode
	Article no.	900416
	Dimensions	Ø 3 mm
	Description	Spatula electrode, flexible
	Article no.	900417
	Dimensions	6 mm 2 mm

EBERLE

EBERLE GmbH & Co. KG

Dachsteinstraße 2
75449 Wurmberg
Germany

Phone +49 (0) 7044 9611-288

sales@eberle-med.com
www.eberle-med.com

